

1

2

Manual of the omnistat

3

Continuous culture system for experimental evolution

4

5



6

7

8

9

written by

10

David Ekkers

11

12

13 Index

14	Components description of the bioreactor system	4
15	Overview of the system.....	4
16	The bioreactors	4
17	The bioreactor	4
18	Module frame.....	4
19	Pumps.....	5
20	Ipcn24.....	5
21	Ismatec Reglo ICC	5
22	Ismatec Reglo digital	6
23	Masterflex Easy Load L/S.....	6
24	Ismatec MS-CA 4/840.....	6
25	Bottles	6
26	Tubing.....	7
27	Influx.....	7
28	Efflux and sampling	7
29	pH control.....	7
30	Chimney.....	8
31	Aeration	8
32	Migration.....	8
33	Antifoam tubing.....	8
34	Filter sterilization tubing	8
35	PH-sensor	9
36	O ₂ -sensor	9
37	Air and N ₂ input	9
38	Stirrers	9
39	Waterbaths.....	10
40	Pinch valve system	10
41	Computer system	10
42	NI CRIO.....	10
43	PC.....	11
44	Security.....	11
45	Software	11
46	Protocols.....	15

47	Making vitamin and metal 1000x stock solutions for CDMPC	15
48	Making filter-sterilized media	15
49	Building up and running the bioreactor	16
50	Handling and calibrating pH sensors	16
51	Handling and calibrating oxygen sensors	17
52	Handling the waterbath	18
53	Setting up gas flowmeters.....	18
54	Setting up the switch system	18
55	Measuring dilution rate.....	20
56	Graphical user interface	20
57	Using Dizanta Compass	21
58	Programming and setting up migration pump.....	21
59	Closing down the bioreactor	23
60	Building up the bioreactor:.....	24
61	After connecting all of the tubing	24
62	Gluing stoppers after autoclaving	25
63		
64		
65		

66 **Components description of the bioreactor system**

67 **Overview of the system**

68 The bioreactor system consists of 24 individually digitally controlled and logged bioreactors. The
69 bioreactors are grouped into 6 modules of four bioreactors that sit together in an autoclavable
70 stainless steel frame. The system is operated via a real time Labview application that is run from a
71 cRIO (NI compact RIO) computer. An emergency power supply in combination with the UPS
72 (uninterruptable power supply) system ensures a continuous power supply to the bioreactor when
73 an experiment is running (Fig 1). The pH in the bioreactors is kept constant with additions of base
74 (NaOH) initiated by pH sensor data from the bioreactors. It is also possible to measure the DO
75 (dissolved oxygen). Bioreactors can be connected via tubes to create corridors between the
76 bioreactor to induce migration between bioreactors (patches). Additionally computer controlled
77 pinch valves provide the possibility to vary the medium components that are added to the
78 bioreactors (Fig. 2).

79

80 **The bioreactors**

81 **The bioreactor**

82 The bioreactors primarily consist of a metal head plate and a 100 ml (standard volume), 150 or 250
83 ml glass reactor vessel (*durant GL45*) (fig. 3). These two components are secured with a hollow
84 screwcap that clamps the glass reactor vessel and the metal head plate with a silicone ring-seal that
85 sits between the two (fig. 4). The head plate has 9 ports (fig. 4): 2 sensor ports (pH, Oxygen or
86 turbation 12 mm \varnothing), 3 static ports (2x 2 mm \varnothing , 1x 3 mm \varnothing) for (chimney, nutrient, nitrogen gas en
87 base additions) and 4 height adjustable ports (efflux, sampling, migration and aeration 2 mm \varnothing) (fig.
88 4). The sensors are secured with a clamping nut and a silicon O-ring (fig. 4). If one of the two sensor-
89 ports is not used one can be closed off with a closing nut and a silicone O-ring (fig. 4). The height
90 adjustable ports consist of a stainless steel tube, that is secured into the port with a PTFE ferrule and
91 a stainless steel clamping nut. If a port is not used the port can be closed off by inserting a solid
92 stainless steel closing rod inside the port instead of the stainless steel tube (fig. 4). For aerobic
93 growth a special aeration rod can also be fitted to the bioreactor to maximize oxygen input in the
94 culture (fig. 4). For homogenization of the culture each bioreactor contains a PTFE coated stirring rod
95 (*VWR 38 mm 442-0398*).

96

97 **Module frame**

98 The bioreactors sit in a stainless-steel frame that holds four individual bioreactors in total (6 modules
99 in total). This frame allows the bioreactors to be autoclaved more easily and secures the bioreactors
100 when they are running. The frame consists of a top-plate and a height adjustable bottom-plate which
101 are both supported by four legs. A bioreactor centralizing disk that sits between the reactor vessel
102 and the screwcap allows for the firm and centralized placement of the bioreactors in the top-plate of
103 the frame (fig. 5+6). The top plate also features cutaways that function as bottle holders (8X), tube
104 guides (16X) and tube-holders for calibration of the sensors (10X). If necessary the bottle holders and
105 tube guides can be closed off to prevent bottles or tubing sliding out of the frame by securing the
106 closing rail with two pins on the module frame (Fig. 5). Beneath the bioreactors the stirring magnet

107 holders are located on the bottom-plate. The bottom-plate is height adjustable with hex-nuts that
108 allow the stirring magnets to be placed at different heights depending on the size of the reactor
109 vessel (100, 150 or 250 ml). A set of tools are available to unscrew nuts for various modifications to
110 the bioreactor system (fig. 7).

111

112 **Pumps**

113 All the influxes and effluxes of liquids that take place inside the bioreactor are controlled via
114 peristaltic pumps. These pumps can all be controlled either manually or via the real-time Labview
115 application to set their speed and direction. More intricate control of the pumps can be programmed
116 with the add-on software (*Dizanta-Compass*). This will be more extensively discussed in the software
117 section of the manual. Each pump uses plastic cassettes to secure the tubes in the pump. These
118 cassettes keep the correct amount of tension on the peristaltic tubing to maximize pumping
119 accuracy and minimize wear on the tubing. The pumps use special 2-stops (IPCN-24) or 3-stops
120 (Reglo ICC and Reglo Digital) tubing. The stops on the tubing serve as guides that allow the tube to be
121 securely fitted in the cassettes. The tubes necessary for the pumps are made of PharMed. This
122 material is resistant to chemicals, autoclaving and physical wear by the peristaltic pumps. The tubes
123 either come with glued non-autoclavable stoppers (PharMed Ismaprene) or welded autoclavable
124 stoppers (PharMed-BPT).

125

126 **Ipcn24**

127 The IPCN-24 is a 24 channel peristaltic pump with a planetary drive with 8 actively driven stainless
128 steel rollers for extremely high accuracy with flowrates between 0,0004-11 ml/min (*Ismatec IPCN-24*)
129 (fig. 8 and 9) . The pump can either operated offline by putting the flow settings in by hand, with an
130 analog signal (not operational) or serial communication via RS232. The IPCN-24 usually controlled by
131 the Labview interface via the RS232-port. In this application the speed and direction of the pump can
132 be set. There are 3 ipcn-24 pumps that each operate 2 modules (8 bioreactors). The pump is
133 primarily used for medium influx as well as efflux. To generate the required dual flow directions, the
134 efflux tube is put in reverse orientation compared to the influx tube. Furthermore the efflux tube
135 always has a bigger inside diameter (usually one or two sizes bigger) than the influx tube to prevent
136 the influx to become greater than the efflux. The influx and efflux take up 16 pump-channels (8x
137 influx + 8x efflux) of the 24 available channels of the pump. Up to 8 additional inputs (1 per
138 bioreactor) can be used for things such as anti-foam, antibiotics, sugars, signal molecules etc. For
139 these additional channels, the pumping speed can only be changed in respect to the dilution rate by
140 the chosen ID of the peristaltic tubing.

141

142 **Ismatec Reglo ICC**

143 The reglo ICC is a 8 roller 4 channel digitally controlled peristaltic pump with independent control of
144 all four channels (*Ismatec REGLO ICC MS-4/8*) (fig. 10 and 11). In the usual setup this pump
145 dynamically adds NaOH solution into each bioreactor to keep the culture pH constant. The pump is
146 connected via a USB to the cRIO and controlled by the Labview application through a PID
147 (proportional-integral-derivative) loop. This loop proportionally activates the pump channels based

148 on the deviation of the measured pH-value in the bioreactor compared to the set-point value. The
149 three parameters (proportion, integral and derivative) can be adjusted in Labview to optimize the pH
150 control. Each module has one 4-channel pump making a total of 6 pumps for the entire system.

151

152 **Ismatec Reglo digital**

153 The reglo digital is a 12 roller 4 channel digitally controlled (RS232) peristaltic pump (*Ismatec REGLO*
154 *Digital MS-4/12*) (fig. 12). It functions as the migration pump continuously transfers culture back and
155 forward between bioreactors at a low rate to create a corridor. To do this the pump alternates the
156 pump direction in time. At what interval the direction is switched and what speed the pump rotates
157 can be programmed with Dizanta Compass. The pump can also be used for dynamical additions of
158 chemicals to the bioreactor such as antifoam, nutrient components, stressors(eg antibiotics,) or
159 signal molecules.

160

161 **Masterflex Easy Load L/S**

162 The Masterflex Easy Load L/S is an offline pump that can be manually set for direction and speed. The
163 pump is used to filter-sterilize media by transferring is through a filter-unit into a sterile bottle. There
164 are two pumps in the bioreactor lab to allow to make two bottles of media in parallel. These pumps
165 do not require any tubing-cassettes or stoppers on the tubing, instead tubing can be directly fed into
166 the pump.

167

168 **Ismatec MS-CA 4/840**

169 Ismatec MS-CA 4/840 Is an offline pump with four channels, 8 rollers that has one pump rate in
170 which only the directional of rotation that can be controlled. The pump can be used in emergency
171 situation where additional pumping is temporally needed. The pump works with the same cassettes
172 and 3-stops tubing as the Reglo pumps.

173

174 **Bottles**

175 All bottles used for the bioreactor system use four-port GL45 screwcaps. Using these caps allows for
176 great flexibility in setting up any experiment. These caps can be screwed on 100 ml up to 20 L bottles.
177 The caps consist of a screwcap with four extruding channels with external thread, four plastic nuts
178 with internal, four silicone ring-seals which sit between screwcap channel and the screw nut to
179 create an air tight seal between the tubing and the bottle. Three stainless steel nuts are secured to
180 the tubing that goes into the bottle to keep it at the bottom. The four ports on the screwcaps are
181 usually used for 3x input or output, 1X air-filter (to either release pressure or vacuum that is built up
182 into the bottle). The bioreactor system uses 20 L medium and waste medium bottles (Fig. 14 and 15).
183 For the medium, brown glass bottles are used to reduce photic reactions in the medium, while the
184 transparent bottles are used for waste medium and storing autoclaved water. The tubing can be set
185 up in a way that one bottle of medium or waste medium supplies either one or a multitude of
186 bioreactors. For other additions such as acid, base or antifoam smaller bottles ($\leq 2L$) are used.

187 **Tubing**

188 The tubing that is used for the bioreactor systems can be divided into two types. Firstly PharMed
189 peristaltic tubing which is used in the peristaltic pumps. This beige tubing is made from a material
190 that allows for prolonged accurate pumping and is very temperature and chemically resistant.
191 Secondly silicone transparent transfer-tubing, this type of tubing is used for all other connections
192 made in the bioreactor system.

193

194 **Influx**

195 The influx tubing starts with a quick-connector that attaches to the medium bottle. From here
196 silicone tubing (ID 2 mm) is attached connecting to the peristaltic tubing (ID 1.30 or 1.02 mm
197 PharMed BPT), the diameter of the peristaltic tubing might vary based on the desired dilution rate of
198 the bioreactor. The other side of the peristaltic tubing is connected to another section of 2 mm
199 silicone tubing that connects to the top end of the liquid-entry-system. The side-input of the liquid-
200 entry-system is connected with a 6 mm silicone tubing to a gas filter which has another piece of 6
201 mm silicone tubing ending in a luer-fitting that will connect to one of the output tubes from the N₂-
202 gasflowmeter. This gas-input functions both as anaerobic ventilation of the bioreactor headspace as
203 well as forcing down the droplets of medium through the liquid-entry-system. The liquid entry-
204 system is a glass drip-feed system for the medium. It has a top-input where the media is fed in, a
205 side-input which connected to a N₂-gas output that pressure forces the media down into the
206 bioreactor via the bottom-output into the medium-port on the bioreactors. The primary function of
207 drip-feed system is to eliminate the chance of the culture to grow back upstream into the medium
208 input when running the bioreactor. This is achieved by dripping in the medium instead of using an
209 uninterrupted medium feed. The gas-input additionally makes sure that the droplets that fall down
210 the drip-feed are pushed into the bioreactor directly. This is important because otherwise medium
211 will accumulate in the bottom of the liquid-entry-system until the weight of the medium overcomes
212 the capillary force within the bottom output, causing a sudden flush of medium in down into the
213 culture causing a more irregular feed into the bioreactor.

214

215 **Efflux and sampling**

216 The efflux tubing starts from the efflux-port on the bioreactor with a short 1 mm silicone tubing that
217 splits up into two tubes, one 1 mm tube ending with a female luer-fitting from where sampling can
218 be done. The other tube changes to a 2 mm silicone tubing that is subsequently connected to a
219 peristaltic tubing (ID 1,52 mm PharMed BPT) that goes into the ipcN-24 pump. This tube is then
220 connected to another piece 2mm silicone tubing that ends in a quick-connect that attaches to the
221 waste medium bottle.

222

223 **pH control**

224 This tubing starts with a connector that attaches to a connector on the base/acid bottle. From this
225 connector 1 mm silicone tubing is used that splits up into four 1 mm silicone tubes. Each of the four
226 tubes is connected to a peristaltic tube (ID 0.89 mm PharMed BPT) which goes into the pH-pump
227 (Reglo ICC). The peristaltic tubing is then again connected to 1 mm silicone tubing that goes all the

228 way to the pH-port on the bioreactor. Great care should be taken when connecting the connector to
229 the base/acid bottle after autoclaving. Clean both connector ends with ethanol and purge the
230 base/acid bottle to remove all the possible air in the tubing. Then attach the tubes sterilely at the
231 flame.

232

233 **Chimney**

234 The chimney is attached to chimney-port on the bioreactor with a 2 mm silicone tubing , this is then
235 converted to a 6 mm tube. This 6 mm tube is filled with a twisted piece of aluminum foil which
236 prevents condensation water in the chimney to be pushed up into the gas filter. From this tube the
237 chimney spits up into two outputs each with a gas filter attached. The one of the two filters serves as
238 a backup filter in case one of the two filters fails.

239

240 **Aeration**

241 The aeration tubing starts with luyer-fitting that is attached to one of the gas outputs of the air
242 flowmeter. From there a 6 mm tube is attached to a gas filter, from the gas filter a 6 mm silicone
243 tube is converted to a 1 mm silicone tube that attaches to aeration-port on the bioreactor.

244

245 **Migration**

246 The migration tubing consists of a 1 mm tube attached to the migration-port on one bioreactor. This
247 tube is then attached to a peristaltic tube (ID 0.89 mm PharMed BPT) which connects to another
248 piece of 1 mm silicone tubing that attaches to the migration-port on the other bioreactor.

249

250 **Antifoam tubing**

251 The antifoam tubing or tubing meant to add any other type of liquid to the bioreactor, starts with a
252 luyer-fitting that attaches to the output of the antifoam bottle from there a peristaltic tube (ID 0.89
253 mm PharMed BPT) is attached. This is then connected again to a 1 mm silicone tubing that attaches
254 to the antifoam-port on the bioreactor.

255

256 **Filter sterilization tubing**

257 This piece of tubing is used to filter-sterilize medium. It does so by submerging one end of the tubing
258 in a bottle of freshly prepared unsterilized medium and pumping (Masterflex easy load) this medium
259 through a filter into a sterile medium bottle. The tubing starts with a long piece of 6 mm silicone
260 tubing that is weight down on the upstream end by two stainless steel bolts that are slipped on the
261 tubing. On the downstream side, the tubing is attached to the input side of a medium filter(). From
262 the output of the medium filter another piece of 6 mm silicone tubing is attached that converts to a 3
263 mm silicone tubing that attached to the input port on the medium bottle.

264

265 **PH-sensor**

266 To measure the pH in the bioreactors during experiments each bioreactor can be equipped with a pH
267 sensor which is connected to a controller. For each bioreactor module (four bioreactors) one
268 controller is used (Fig. 15 and 16). The controller receives the signals from each pH-sensor which is
269 then displayed real-time on the screen of the controller and subsequently send using 0-20 mA signal
270 to the bioreactor computer. This pH value is then logged and used to control the pH pump via a PID-
271 loop.

272 The pH-sensor is a solid gel pressurized sensor with a silver/silver chloride reference suitable
273 autoclaving (*Applisens pH sensor for mini bioreactor 8 mm 15 cm*) (Fig. 17). The sensor typically
274 withstands around 20 autoclaving cycles before it should be replaced. After autoclaving and
275 calibration a sensor will remain accurate for up to 3-4 weeks if it cannot be recalibrated during its
276 use.

277

278 **O₂-sensor**

279 To measure the DO (dissolved oxygen) in the bioreactors during experiments each bioreactor can be
280 equipped with a DO-sensor which is connected to a controller (Fig. 18). The controller receives the
281 signal from the DO-sensor which is then displayed real-time on the screen as either a percentage or
282 amount (mg/l). The sensor signal is subsequently send using 0-20 mA signal to the bioreactor
283 computer. This value is then logged and can also be used for controlling other systems such as pinch
284 valves to change gas mixtures.

285 The DO sensors used for the bioreactors are polarographic (*Applisens 8 mm 15 cm polarographic DO₂*
286 *sensor for mini bioreactor*) (Fig. 19). The sensors measure partial pressure of dissolved oxygen in
287 liquid. The sensor is equipped with a gas-permeable PTFE membrane that allows oxygen to pass to
288 be measured by the sensor. The sensor operates by the Clark-cell principle: a fixed polarization
289 voltage of 675mV across the sensors anode/cathode reduces O₂ molecules. This reduction results in a
290 current that is linearly dependent on the amount of reduced O₂ molecules. The sensor has a
291 protective cap covering the membrane to avoid mechanical damage and to limit evaporation of
292 electrolyte.

293

294 **Air and N₂ input**

295 The air and N₂ input system consists of a (1) N₂-supply and an air-supply with both a reduction valve
296 with a bar meter (Fig. 21 and 23), (2) distribution tubing from the reduction valve (Fig. 20) and (3) 4-
297 channel flowmeters with a main open/close valve and accurately adjustable needle valves for each
298 specific channel to regulate gas input into the bioreactors (Fig. 22 and 24). The range for the air-
299 flowmeters is 0-100 ml/min and for the N₂-flowmeters 0-25 ml/min.

300

301 **Stirrers**

302 The bioreactor system uses magnetic submersible stirrers (*Cimarec micro*) that are placed inside of
303 the module frame stirrer holders (fig. 26) after autoclaving the bioreactors. The stirrers are then

304 submerged into the water bath together with the bioreactors. The magnetic stirrers propel stirring
305 rods inside the bioreactors to agitate and homogenize the culture. The magnets are offline and
306 operated per four stirrers via a hand operated tuning-knob located above the bioreactor (fig 25).

307

308 **Waterbaths**

309 The bioreactor are being kept at a constant temperature by submerging the bioreactor module frame
310 into a temperature controlled waterbath (*Julabo MB*) (Fig. 27). The waterbaths are hooked up to the
311 bioreactor computer via a RS232 connection and can be either controlled via the bioreactor or
312 manually.

313

314 **Pinch valve system**

315 The pinch valve system is designed to be able to change liquid or gas input or output in time within
316 the bioreactors. This is achieved by an array of computer controlled pinch valves that can open or
317 close tube channels connected to the bioreactor. Examples for the use of such a system can be
318 varying aerobic versus anaerobic conditions in time by changing gas input in time. Adding stress
319 molecules to the bioreactor to create evolutionary bottlenecks in time. Changing media components
320 during a run to induce metabolic switching. Because the pinch valve system is hooked up to the
321 labview application pinch valve activation can also be made dependent on the output parameters of
322 the bacterial culture such as acidification rate or turbidity. The pinch valve system consists of three
323 switch boxes (one for every two modules) each containing an array of 8 double channel pinch valves
324 (*Sirai 3/2 NC-NO solenoid pinch valve*) (Fig. 30). The pinch valves on the switch boxes can be
325 controlled manually or from the cRIO. From the cRIO the pinch valves can either be operated directly
326 from labview application or indirectly via the add-on software Dizanta Compass. The latter one
327 allows the pinch valves to programmed with more complex logic behind them. The control signals
328 from the cRIO are send from an i/o card (Scheme 1) that it translates into 0-5 V TTL signals that go to
329 three control boxes located on the wiring guide above each two bioreactor modules (Fig. 28 and 29).
330 The signals are then amplified and transmitted to the pinch valves. In an 'deactivated state' a spring
331 in each pinch valve always keeps one of the channels closed while the other channel is open. When
332 the valve is activated, the coil inside the valve opens up the formerly closed channel and at the same
333 time closes off the other formerly open channel.

334

335 **Computer system**

336 **NI CRIO**

337 The operating system of the bioreactor is controlled by a real time target, a cRIO (*National*
338 *instruments compact RIO-9067, 8 slot integrated dual-core controller and artix-7 FPGA*) computer
339 that is mounted via a DIN-rail to the wall. The computer is powered by a 24V DC output converter (*NI*
340 *PS-15 power supply*) that is connected to the standard 220V AC (Fig. 31 and 32). Sensor data from the
341 pH and DO (dissolved oxygen) sensors is converted by the sensor controllers into a standard 0-20 mA
342 signal. These signals enter the computer via two current input modules (*NI 9208 24-bit current input*
343 *module with D-sub*). Pinch valves are operated from the computer via TTL digital output module (*NI*

344 9403 with D-sub 32 channel TTL digital input/output module). The pumps controlling the pH (REGLO
345 ICC) inside the bioreactors are controlled via a serial signal send from the USB port from the cRIO
346 that is branched into 6 channels via a USB hub located on the wiring guide above the first bioreactor
347 module (Fig. 33). The migration pumps(REGLO digital), flux pump (IPCN-24) and the water baths are
348 connected to the computer via three serial RS232 modules with 4 channels each (NI 9870 4-port
349 RS232 serial module W/4 10P10C-DE9). For an overview of all connections see the wiring scheme
350 (Scheme 1).

351

352 PC

353 The cRIO is linked to a PC. During the running of the bioreactor system the software is always
354 operated by the cRIO which holds a database with the most recent data. The cRIO updates the PC
355 four times every second with the latest data by messages. The PC updates the tags in the database to
356 the controls available on the GUI (graphical user interface). The cRIO data that is logged is
357 transferred to the desktop computer (D:ProControl/data/).

358

359 Security

360 The whole computer system is connected to an emergency power supply that will kick in a matter of
361 1-10 sec when there is a power failure. Furthermore all the other systems except for water bath 4, 5
362 and 6 are connected to and uninterruptible power supply (UPS) that can buffer a power dip of up to
363 40 sec long. If for some reason the desktop computer is turned off or crashes while the bioreactor
364 system is operating the cRIO automatically logs the data on to an external hard disk that is connected
365 to the usb port via a hub (Fig. 33). When the desktop computer is operational again the data will be
366 saved on the computer again.

367

368 Software

369 The bioreactor system uses a Labview based GUI (graphical user interface). The first tab of this
370 interface allows the user to set the input values of the pumps (activity, speed and direction),
371 waterbaths(temperature), pH (set point and PID) and pinch valves (Fig. 34 and 35). The second tab of
372 the GUI plots all logged parameters such as, pH, DO, temperature, pinch valve activity, pump speed,
373 pump directionality and turbidity (not operational) per module or bioreactor (Fig. 36). This allows for
374 monitoring of the bioreactors during running of an experiment as well as adjustment of its
375 parameters. Although the Labview GUI allows easy reliable bioreactor use it does not allow the user
376 to easily program more complex logic between or on top of the input and output parameters.

377 Software framework

378 The software framework is a database subscription framework. The compact RIO (cRIO) from
379 National Instruments (NI), holds the database with the most recent data. The cRIO updates the PC
380 with the latest data four times per second by messages. The PC updates the tags in the database to
381 the controls available on the graphical user interface (GUI). The tag configuration sets which tags are
382 available in the system for data acquisition and data generation. Tags are clusters, which hold
383 information about the connection type, ranges, alarm ranges, scaling etcetera.

384 To minimize system failure the compact RIO runs independently from the PC application. Therefore
385 the cRIO application is running while powered. When a working LAN connection has been
386 established, data will be transferred and the GUI will be updated. By default the received databases
387 are stored in the format TDMS (Technical Data Management Streaming) on the host PC. The PC
388 application can run and function while the cRIO is not connected or turned off, though updates are
389 delayed until the cRIO is reconnected or switched on.

390 The GUI shows the latest version of the database to the user. There are indicators to visualize data
391 and controls which both visualizes the data and can change values. The controls (and indicators) can
392 be of the type Boolean which allows for switching between TRUE and FALSE, which represents for
393 example an ON and OFF state or OPEN and CLOSE state. Each control represents a tag from the
394 database and the value is updated on the GUI. The tag contains properties described in the tag
395 configuration. This can e.g. limit a control to a range of input values. In order to access the other tag
396 properties, the context menu shows the information stored. By selecting “View parameters Tag”, the
397 popup windows allows for visualization of the tag properties and editing of writable tag properties.
398 The same holds true for editing and viewing the properties of the PID settings. When editing,
399 pressing ‘Cancel’ will revert the changes and ‘Ok’ will store the new parameters.

400

401 PID loops

402 Any set of parameters can be turned into a PID loop. A PID loop has a setpoint (SP), process variable
403 (PV), output parameter (OP), control switch (CTR), output rate, output range and a default value
404 when switched off. The P stands for proportional gain, I stands for the integral time in minutes and
405 the D stands for the derivative time in minutes. The concept is that the system changes the output
406 action in order to minimize the error between the set point and the process variable. The error at
407 time point t , $e(t)$, is calculated from

408

$$409 \quad e(t) = SP(t) - PV(t) \quad \text{eq. 1}$$

410

411 The controller action at timepoint t is as following

412

$$413 \quad u(t) = K_c (u_p(t) + u_i(t) + u_d(t)) \quad \text{eq. 2}$$

414

415 where K_c stands for the controller gain (P). The proportional action u_p is calculated by

416

$$417 \quad u_p(t) = e(t) \quad \text{eq. 4}$$

418

419 The integral action is calculated by

420

$$421 \quad u_i(t) = 1/T_i \int_0^t e dt \quad \text{eq. 4}$$

422

423 where T_i stands for the integral time (I) in minutes. The derivative action is calculated by

424

$$425 \quad u_d(t) = T_d \frac{de}{dt} \quad \text{eq. 5}$$

426

427 where T_d stands for the derivative time (D) in minutes.

428

429 Configuration

430 The configuration tool allows the administrator to modify the default tag values and default PID
431 settings. For the tags:

432 Tag name: tag name

433 Units: units of tag

434 IO Type: type of connection

435 Access: read only (R), write only (W), read and write (RW)

436 Alarm LL, L, H, HH: alarms values for low low, low, high and high high alarms

437 Auto ack: auto acknowledgeable

438 Raw: range of raw values

439 Eng: scaling of raw values to engineering values

440 Connection: connection address

441 Description: extra information for user or device

442 Design: hard limits for range of setpoints

443 Default: start up value

444 For the PID settings:

445 P: proportional gain

446 I: integral time in minutes

447 D: derivative time in minutes
448 Output rate: rate for the output in engineering units per minute
449 Range output: range of the output
450 Default value: default value when switched off
451 Setpoint: setpoint
452 Process variable: process variable tag name
453 Output parameter: output parameter tag name
454 Control: control tag name

455

456 Dizanta Compass

457 For more complex programming the add-on software Dizanta Compass is used. Dizanta Compass
458 allows you to take any of the input or output parameters from the bioreactor system and program
459 them into a logic scheme that produces input values for the GUI of bioreactor system (Fig. 37). When
460 the program is run, it simply acts as a virtual user changing input values in the labview GUI according
461 to what is programmed in Dizanta. This software allows for a very flexible and variable use of the
462 bioreactor system. An example of logic on top of the parameters could be a preprogrammed
463 temporal variation of chemical additions to the bioreactor, by preprogramming the pinch valves and
464 pumps to add certain type and amount of chemicals through time to the bioreactor. An example of
465 logic between parameters could be making the dilution rate of the bioreactor dependent on the
466 growth speed of the culture, by making the bioreactor flux-rate dependent on acidification rate or
467 turbidity of the culture.

468

469 Protocols

470 Making vitamin and metal 1000x stock solutions for CDMPC

- 471 1. Start by adding 50% of the end volume of demi water to a flask and adding a clean
472 stirring rod.
- 473 2. Next weight and add all vitamins or metals to the solution while it is stirring. See scheme
474 2 for all components.
- 475 3. Then add water until the 100% volume is reached.
- 476 4. For the vitamin solution the pH must be adjusted in order to dissolve all vitamins. Add
477 NaOH solution until a pH of 7 is reached. Now the components will dissolve in a matter
478 of minutes. As the vitamin solution dissolves it will acidify the solution, so the pH must be
479 adjusted several times until all vitamins are dissolved.
- 480 5. Then the 1000x vitamin or metal solution is filter sterilized by passing it through a 0,2 μm
481 filter. Aliquot the solution into 20 ml portions for making 20 L medium at once and
482 aliquot smaller volumes for making smaller amount of medium. The aliquots stock
483 solution should then be frozen at -20 C. The stock solutions should only be defrosted
484 once when it is used!

486 Making filter-sterilized media

- 487 1. To make for a filter-sterilized medium such as CDMPC (chemically defined medium for
488 prolonged cultivation) you start with a 20 L medium bottle and a 20 L transparent medium-
489 preparation bottle.
- 490 2. The medium bottle should have already the filter-sterilization tubing attached to the input
491 channel before autoclaving. This tubing is roiled up into a packed of aluminum foil and taped
492 against the side of the medium bottle to make it easy to autoclave the whole thing. This
493 setup makes sure that the tubing downstream of the medium filter is completely sterile
494 without risk of contamination because there is no need for connecting any tubes. Make sure
495 that the medium bottle is height marked at the 20 L mark with a piece of tape.
- 496 3. The medium-preparation bottle should be filled up with ± 18 L of demi water and a big
497 stirring rod should be added to the bottle. The water is autoclaved in one of the big
498 autoclaves on the fourth floor.
- 499 4. Now weight the medium components (see scheme for CDMPC recipe). Add them to the
500 medium-preparation bottle with 18 L water while the liquid is stirred with a stirring rod. The
501 components dissolve fasted when the liquid is stirred at the highest possible speed.
- 502 5. In the case of making CDMPC the amino acids should be added first, wait until everything is
503 dissolved, this can take up to 10-45 min depending on the temperature of the water.
- 504 6. Next you add the buffers and carbon sources, again wait until everything is dissolved.
- 505 7. The vitamins and metals can be added from a 1000X or 500x frozen stock solution.
- 506 8. When all components are dissolved in the medium, the pH should be adjusted (pH 6.5 in the
507 case of CDMPC). To do so the pH sensor is hanged into the medium while stirring at a low
508 speed (make sure to regularly calibrate the sensor!). Then NaOH or HCl solution is added
509 until the pH is stably set at the desired value.

- 510 9. Now the media is ready to be sterilized. The medium-preparation bottle is put on the bench
511 next to the Masterflex pump. The empty and sterile medium bottle is put on the ground
512 below the pump.
- 513 10. Now the packet containing the filter-sterilization tubing is opened and the upstream part of
514 the tubing is lowered into the medium preparation bottle. Use a piece of tape to stick the
515 auto-going tubing to the side of the medium-preparation bottle, this will prevent the tube
516 from falling out of the bottle during pumping. From the medium preparation bottle the
517 tubing is led through the Masterflex pump and clamed into the rollers. Make sure that the
518 medium filter is always hanging in an upright position under the Masterflex pump. This is
519 important to effectively purge any trapped air from the filter via the degassing valve.
- 520 11. Before activating the pump check if the tubing is set up correctly, check for any places where
521 the tubing might be folded. Then open up the degassing valve on the medium filter and
522 activate the pump.
- 523 12. During the first minute of pumping make sure that all the air is purged from the filter. Air
524 tends to get stuck in between the lamella of the filter, firmly tap the filter to release these air
525 bubbles.
- 526 13. When almost all of the ± 18 L of medium is pumped into the medium bottle more water
527 should be added. Make sure to stop the pump before air is sucked into the tubing. Then add
528 1L of demiwater and swirl the water around the bottle to make sure a much of the medium
529 components are added to the medium. After the 1 L is pumped through add another 1 L of
530 demi water until the 20 L mark on the medium bottle is reached. Now clamp off the 3 mm
531 silicone tubing, that is attached to the medium bottle input with two clamps. Cut the tubing
532 just upstream of the two clamps to de-attach the medium filter from the medium bottle.
533 Write down the preparation date on a tape on the medium bottle. The medium bottle is now
534 ready to use. CDMPC stays good for at least a month.
- 535 14. Make sure to immediately clean the medium preparation bottle and the filter-sterilization
536 tubing with NaOH solution and demiwater directly after use to prevent any kind of microbial
537 growth!!! After cleaning, autoclave a new batch of water in the medium-preparation bottles.
538 These can then be stored ready to be used for preparing the next batch of media. The filter-
539 sterilization tubing should be straight away fitted with a new medium filter, repackaged into
540 aluminum foil and autoclaved with a cleaned medium bottle. If this is not possible the pieces
541 of tubing should be hanged to dry.
- 542

543 **Building up and running the bioreactor**

544

545 **Handling and calibrating pH sensors**

546 The sensors should be calibrated before autoclaving the bioreactors in buffer solutions of pH-7
547 (green) and pH-4 (red) .

- 548 1. The tubes with calibration solution should be placed inside the bioreactor frame so that they
549 are submerged into the waterbath. Make sure the waterbath has the same temperature as
550 the running temperature of the experiment (normally 30 °C).
- 551 2. Remove the pH sensor cap and clean the sensor tip with demi water (always handle the
552 sensors with gloves!).

- 553 3. Put the sensor inside of the red 4 pH buffer solution.
- 554 4. Wait for 15 min for the buffer solution and sensor to attain the same temperature as the
- 555 waterbath. Then rinse the sensor with demi water, make sure you keep the sensor in the
- 556 upright position to prevent movement of sensor gel upward in the sensor and bubbles from
- 557 forming inside the sensor tip.
- 558 5. Enter the menu of the sensor controller by pressing the menu key from the main
- 559 measurement screen and select 'calibration'.
- 560 6. Select the pH channel that you want to calibrate and press enter.
- 561 7. Enter again and check if the pH value of the sensor is stable and adjust the value to the pH of
- 562 the buffer solution (pH4 for the red solution) if needed.
- 563 8. Then transfer the sensor to the second buffer solution, make sure to first rinse the sensor
- 564 with demi water before putting it in the second solution. Again press enter when the signal is
- 565 stable and adjust the value to the buffer value (pH 7 for the green solution).
- 566 9. Slope and offset value are now displayed and saved in the calibration log.

567 A calibrated sensor remains usable for 3-4 weeks depending on the required accuracy. After use in an
568 experiment the sensor is again autoclaved with the bioreactor and subsequently removed from the
569 sensor ports and rinsed thoroughly and stored in sensor solution until the next use. Sensors can be
570 autoclaved around 20 times before they tend to fail. Sensor failure will be indicated by the failure to
571 calibrate.

572

573 **Handling and calibrating oxygen sensors**

574 The sensor can be best calibrated after sterilization because the autoclaving can alter the sensitivity
575 of the sensor. After sterilization the sensor should at least be polarized for 6 hours before calibration.
576 The DO sensors are mounted inside the bioreactor with a sensor clamping nut.

- 577 1. Remove the O₂ sensor cap and clean the sensor tip with demi water (always handle the
- 578 sensor with gloves!).
- 579 2. Visually inspect the sensor's PTFE membrane. It should be undamaged and clean.
- 580 3. Carefully slide the sensor through the sensor port of the bioreactor. Make sure that the
- 581 clamping nut is sufficiently loose but still attached to the bioreactor cap and make sure that
- 582 there is an undamaged silicone O-ring at the sensor port.
- 583 4. Gently lower the sensor in its correct position.
- 584 5. Clamp the nut tight with a nr. 11 wrench until the sensor is secure, don't tighten too strongly!
- 585 6. Make sure that the protective cap on the top is properly screwed on to the sensor before
- 586 autoclaving.
- 587 7. When the autoclaved bioreactor is back into the waterbath remove the protective cap and
- 588 connect the sensor cable to the sensor, store the protective cap in the designated container.
- 589 Now the sensor should be polarized for at least 6 hours.
- 590 8. After six hours of polarization calibration should be done inside the sterile demiwater in the
- 591 bioreactor. Make sure that the bioreactor is already stabilized at process conditions (stirring
- 592 speed, temperature and gas flow).
- 593 9. Calibration can either be done via a one point (100%) or two points (100% and 0%). If you do
- 594 not need very high accuracy, for instance when you grow aerobic, a one point calibration will

595 do fine. However when you grow anaerobically or semi-anaerobically it is a good idea to do a
596 two-point calibration. For the two point calibration check the manual of the DO control box.
597 10. One-point calibration: saturate the bioreactor at process conditions for approx. 30-40 min
598 with air (100 ml/min). When a stable DO signal is reached perform the 100% calibration on
599 the DO controller by pressing the 'sens' button for 3 seconds. Now the message 'sens. Cal'
600 appears followed by the actual value. Now changes the value of the sample concentration to
601 100% with the 'sens' and 'set1' button. Confirm and end the calibration by pressing the 'set2
602 ent' button. If everything has gone alright the message 'update' will appear.

603

604 **Handling the waterbath**

605 For a good operation of the waterbaths they should be filled with 7 L of demi water(always use
606 demi-water to prevent formation of calciumcarbonate!). Then disinfectant (*aqua resist*) should be
607 added (1.5 ml/L) to prevent microbial growth inside of the waterbath. This antimicrobial solution
608 contains a blue indicator colorant. New antimicrobial solution must be added to the waterbath water
609 whenever the blue color is almost invisible in the water. Make sure to wear gloves when handling the
610 antimicrobial solution or the waterbath water, since the solution is toxic! Because the heaters of the
611 waterbath use a lot of power, care should be taken not to overload the fuse box. When running more
612 than 3 modules make sure not to activate all the waterbaths at once, but instead turn them on three
613 or less at a time, wait for them to reach their target temperature and turn on the remaining
614 waterbaths. When running the system the waterbath need to be topped up with water every other
615 day to prevent the water level to go below minimum. However it is advised to top up the waterbaths
616 every day as good practice. This will prevent problems if one day you happen to forget to top up the
617 water.

618

619 **Setting up gas flowmeters**

620 The gas systems are fed from the air-supply and N₂-supply (5-6 bar) that can be opened or closed
621 with a valve located next to the door of the lab. These valves should normally always be open. From
622 here the gas can be tapped off from the reduction valves located underneath the shelf midway on
623 the bioreactor lab table (Fig. 21 and 23). The gas pressure should be set at 1 bar for a proper
624 operation of the flowmeters. This can be achieved by first opening all open/close valves of the
625 flowmeters that you intent to use during the experiment. Next you gently open the general gas valve
626 located on the side of the reduction valve (Fig. 21 and 23). If you see that the pressure meter already
627 goes up beyond 2 bar while opening the valve, you have to twist the big nob to reduce the pressure
628 because too much pressure can damage the flowmeters. You can now read the current pressure that
629 comes out of the reduction valve into the system. Then adjust the barrage of the output to 1 bar by
630 gently twisting the big nob until pressure meter shows exactly one bar. Finally the amount of added
631 gas can then be adjusted on the channel valves on the flowmeter (Fig. 22 and 24).

632

633 **Setting up the switch system**

634 The pinch valves can be operated manually or via the computer. Each pinch valve has its own
635 tumbler switch that can be put is three positions: the 'middle position' means that the pinch valves

636 are operated via the bioreactor software, switch in the 'left position' is a manual inactivation of the
637 pinch valve (cannot be overwritten by the computer) and switch in the 'right position' will manually
638 energize the pinch valve. Additionally an orange indicator light is located next to the switch indicating
639 whether the valve is energized. To use the pinch valves the control boxes should be turned on by
640 pushing the green switch on the box (Fig. 28). Make sure that all the switches on the switchbox are
641 put into the 'middle' position (software controlled). The tumbler switches on the switch box have a
642 safety lock to prevent unwanted movement of the switches when running. To move the tumbler
643 switch the knob should first be lifted up to unlock. When a pinch valve is activated (manually or
644 software) it first receives a high energetic pulse (activation current) to make the switch movement,
645 after this initial switch a lower current (hold current) will keep running through the activated switch
646 to keep the valves in its position for as long as the switch is on. This feature will prevent the pinch
647 valves to overheat when they are activated for long periods of time. The amount of 'hold current'
648 that is run through the pinch valve can be set with a tuning knob on the control box (Fig. 28). Each
649 pinch valve has two channels: a passively open channel that is closed when the pinch valve is
650 activated and a passively closed channel that is open when the pinch valve is activated. Whenever
651 there is asymmetry in the switch regime you want to induce, always put the more frequently used
652 input in the passively open channel to minimize activation time and thus overheating of the pinch
653 valves. For example, if you would add solution-A for 30 minutes and solution-B for 5 minutes, the
654 input of solution-A should be in the passively open channel of the pinch valve and the input of
655 solution-B should be in the passively closed channel of the pinch valve.

656 A special type of silicone tubing that has relatively soft walls should be used for the pinch valves:
657 HelixMark silicone platinum treated tubing with a 1,98 ID and 3,18 mm OD with a Shore-strength of 50
658 A (VWR 228-1067) or a tubing with similar properties. Never use a type of tubing in the pinch valves
659 which does not have these specifications, because the tube might not close completely or the pinch
660 valves might be damaged! Since this type of tubing is not very rigid it is not recommended to use
661 more than a little piece in the place where the tubing runs through the pinch valves. Using this tubing
662 for longer parts of the tubing might run the risk of tubing folding in on itself and blocking the flow.

663 Depending if you want to switch between gases or liquids, different pinch channel setups are
664 required. If you want to switch between two liquids, you put each tube corresponding to the
665 respective liquid in one of the two channels of the pinch valve. The valve will always open up one
666 channel and close-off the other, enabling you to switch between two liquids with one pinch valve.
667 However if you switch between two gases you will need two pinch valves. This is because gas is
668 always flowing from each channel regardless whether it is closed or open. This means that if a
669 channel to the bioreactor is closed, a venting channel should be opened allowing the unused gas to
670 escape until it is connected to the bioreactor again. This is achieved by splitting up each gas input
671 into two channels that each are fitted inside a channel of the pinch valve. One channel exiting the
672 pinch valve goes to the bioreactor (input channel) and one channel exiting the pinch valve is open to
673 the air (the venting channel). In conclusion whenever you are running the bioreactor, one gas will be
674 connected to their input channel, while the other gas will be connected to its venting channel. This is
675 reversed for both pinch valves when the gas input is switched to the other gas. So two pinch valves
676 are needed to switch between two gases.

677 **Measuring dilution rate**

678 If flow measurements are desired during the running of the bioreactor a Y-section can be attached
679 just upstream of the peristaltic tubing. A piece of 50 cm 2 mm silicone tubing is attached to the Y-
680 section connection and clamped off at the other end. After autoclaving, the bottom of a sterile 2 ml
681 disposable pipet should be sterilely connected to the clamped of piece of silicone tubing. This is best
682 done at the flame by cutting the tube with a sterile scissor. Make sure to additionally clamp off the
683 tube just after the Y-junction. Additionally a 2 μ m syringe filter should be attached to the top end of
684 the pipet (also sterilely!). To measure the dilution rate of a bioreactor the pipet should be secured in
685 a vertical position with one of the rod holding systems. Then the tube attached to the pipet should
686 be opened to the influx channel by removing the clamp. Subsequently medium should be drawn all
687 the way up the pipet to the level just below the little cotton air filter inside the pipet (make sure not
688 to make the cotton wet because it will render the pipet useless for measuring!). Then the medium
689 input attached to the medium bottle close to the y-section should be clamped off. As a result the
690 medium pump will start pump just the medium that is in the pipet. The dilution rate can now be
691 measured by measuring how much time it takes before the pump has pumped away a certain
692 amount (usually timing 1.5 ml is easiest) of medium from the pipet. When finished, the tube
693 attached to the pipet should be clamped off again and the clamp on the media input should be
694 removed.

695

696 **Graphical user interface**

697 To start the operating software on the bioreactor click on the labview icon "PC281 RUG Labview"
698 which is located on the desktop of the PC. You will now see the 'control tab' of labview GUI with all
699 the set point controls visible and ready to use. If somehow the system is not responding as expected
700 you can open NI MAX in order to check the status of cRIO. Within NI MAX you open "remote
701 systems" and choose "NI cRIO 9067-030703de" (Fig. 40). If all seems well in NI MAX, restarting the
702 cRIO often solves the problem. You can do this by pressing the restart button on the top of the
703 screen, a restart takes approximately 10 min to complete.

704 To set the pH for the bioreactors you type in the desired set point and click on the 'play-button' (Fig.
705 35). Make sure that the directionality of the pump is set correctly in respect to the tubing setup.
706 Clicking on the arrow will change the pump direction. The default direction is set at clockwise. When
707 the play button is pressed the PID loop will start to activated the pH pump based on the deviation of
708 the measured pH from the set point. The resulting activity is based on the values that are given to P, I
709 and D, which can be retuned whenever a new type of experiment is started with different
710 acidification or alkalization characteristics. To set the PID value you should right click on the pH
711 window and choose "View parameter Tag". In some cases it might be preferred to use the pH pump
712 for some other function. In this case the play button should not be pressed but the a speed value (%)
713 should be typed in directly in the GUI.

714 Both the medium pump (ipcn-24) and the migration pump (Reglo digital) can only set for speed and
715 direction from the Labview application. When the migration pump is used to create connectivity
716 between bioreactors this is programmed by pumping culture back and forward between two
717 bioreactors. This is achieved by programming the pump to change direction at a set interval using
718 Dizanta Compass.

719 From the GUI the pinch valves can only activated by clicking on them and again clicking to deactivate.
720 However in most cases pinch valves will be used either in a preprogrammed temporal scheme or be
721 made dependent on other variables from the bioreactor. To program such behavior Dizanta Compass
722 will be needed which will be further explained in the following section.

723

724 **Using Dizanta Compass**

725 Dizanta Compass is a licensed add-on software kit that allows for complex logic to be programmed
726 on top of, and between the input and output parameters from the bioreactor. Dizanta compass has a
727 limit of 350 program lines within the current license. To start the program first you need to log in
728 Dizanta Gears. To make a new program you open the editing screen. On the right bottom you can
729 see the amount of available lines. Any program is build up out a number of steps that contain a
730 certain number of actions (Fig. 37). A program is made by selecting one of the 26 preprogrammed
731 input elements (Fig. 38). The input elements are either stand-alone (step, wait, repeat, jump to step,
732 etc.) or relate to input and output parameters from the GUI (activate/deactivate, increase/decrease,
733 set point, not, and, or, etc.) (Fig. 39). When the type of input element is selected the relevant
734 characteristics can be edited (name, input parameter, set point value, time etc.) on the left side and
735 subsequently added to the program (Fig. 38 and 39). After you finished a program for the bioreactor
736 you should test run it within the program to check for bugs. Be aware that sometimes steps with a lot
737 of input elements should be repeated a few times or have additional wait steps of a few seconds to
738 avoid the program to skip certain actions because it runs to quickly through them. If all seems fine
739 you can save the program and open it in the “control” window to run it.

740

741 **Programming and setting up migration pump**

742 To run an experiment with spatial structure, the migration pumps are used to provide connectivity
743 between spatial patches (the bioreactors). The rate of migration within the bioreactor system is
744 expressed as % flux compared to the bioreactor flux. For example a bioreactor with a flux (dilution)
745 rate of 500 $\mu\text{l}/\text{min}$ with a 5% migration rate would transfer 25 $\mu\text{l}/\text{min}$ of culture to the other
746 bioreactor and vice versa in the case bidirectional migration. In order to induce any migration rate in
747 the bioreactor system two values should be calculated: the migration pump rate (R_{mp}) and the
748 migration switch time (T_s) (the time interval at which the pump changes the pump direction). These
749 values depend on four parameters, corridor tube volume (V_c), bioreactor dilution rate (R_d), migration
750 rate (R_m) and switch interval based on growth rate percentage (P_d). To calculate these values see the
751 formula below this section.

752 The volume of the tubing that comprises the corridor between bioreactors should be as small as
753 possible. This is important because selection within this tubing is expected to be different from those
754 in the patches (bioreactors), so you want the smallest possible fraction of cells experiencing these
755 ‘corridor conditions’ compared to the bioreactor conditions. 0,89 mm peristaltic tubing is used
756 because it is the smallest possible PharMed BPT tubing available which is still compatible with the
757 smallest connectors. 1 mm silicone tubing is used for the rest of the tubing length. When making
758 migration tubing it is important to keep the tubing as short as possible to reduce the volume and to
759 make each ‘corridor’ tubing exactly the same so that experiments will induce all the same conditions.
760 When this tubing is finished it is important to accurately calculate the total volume of corridor tubing

761 using the formula stated below. This is important because you need to know the amount of 'dead
 762 volume' that must be pumped every time the direction of the pump is switched before effective
 763 migration is possible. This dead volume is added to the desired 'netto' migration volume to calculate
 764 the total migration pump rate.

765 To achieve bidirectional migration, the migration pump needs to alternate the pump direction at a
 766 certain interval. The switch interval is optimized at 50% of the growth rate in the bioreactors,
 767 meaning that within one generation the pump must have pumped two intervals in each direction.
 768 This value is derived from a tradeoff between switching as little as possible to reduce the impact of
 769 the 'dead volume' on the total pumping activity and on the other hand switching often enough for
 770 cells in each patch (bioreactor) to experience migrated cells within each generation. Based on the
 771 equations both the switch time as well as the migration pump rate can be calculated. It is important
 772 to note that if the flux rate of the bioreactor is adjusted during an experiment both the pump rate as
 773 well as the switch rate should be recalculated in order to ensure a constant representative migration
 774 rate.

V_{pt}	= Peristaltic tubing (μ l)
V_{st}	= Silicone tubing (μ l)
V_{mr}	= Migration rod (μ l)
V_c	= Corridor volume (μ l)
L_{pt}	= Peristaltic tubing length (mm)
L_{st}	= Silicone tubing length (mm)
L_{mr}	= Migration rod length (mm)
D_{pt}	= Peristaltic tubing diameter (mm)
D_{st}	= Silicone tubing diameter (mm)
D_{mr}	= Migration rod diameter (mm)
T_s	= Switch time (min)
R_d	= Dilution rate (μ l/min)
R_m	= Migration rate (μ l/min)
R_{mp}	= Migration pump rate (μ l/min)
P_m	= Migration percentage (%)
P_d	= Switch interval based on growth rate percentage (%)

$$V_{pt} = L_{pt} \times \left(\left(\frac{D_{pt}}{2} \right)^2 \times \pi \right)$$

$$V_{st} = L_{st} \times \left(\left(\frac{D_{st}}{2} \right)^2 \times \pi \right)$$

$$V_{mr} = L_{mr} \times \left(\left(\frac{D_{mr}}{2} \right)^2 \times \pi \right)$$

$$V_c = V_{pt} + V_{st} + V_{mr}$$

$$R_m = R_d \times P_m$$

$$T_s = 60 / \left((R_d / 1000) \times P_d \times \frac{1}{2} \right)$$

$$R_{mp} = (R_m \times 2) + (V_c / T_s)$$

Closing down the bioreactor

At the end of an experiment the module with bioreactors is first disconnected from all in and outputs and then autoclaved in same way as at the start of the experiment. The autoclave module should then be first thoroughly cleaned taking apart the entire bioreactor. Peristaltic tubing that can potentially be used another time should be soaked in 1% NaOH solution overnight and thoroughly rinsed with demiwater. After six weeks all peristaltic tubing should be discarded except for the pH pump tubing this can be used for approx. 18 weeks before replacing. Any dirty silicone tubing should also be replaced. Larger parts of the bioreactor such bioreactor, head-plate, module frame, screwcap etc, can be cleaned in one of the washing machines on the fourth floor. After thoroughly cleaning all parts of the bioreactor they should be immediately dried with pressurized air and then be put in the drying oven on the fourth floor. The bioreactor and modules should then be reassembled and put back in their position ready for use.

Building up the bioreactor:

1. Before building the bioreactor it should be clear how many inputs and outputs are needed and which sensors are required for the specific application. The most important question for the setup of the bioreactors is if you want to run anaerobically or aerobically. If you want to run aerobically extra functions should be added: an air-input, a DO sensor (optional) and possibly an extra input for additions of antifoam solution to prevent clogging up of the chimney.
2. When the setup is clear, add all necessary rods into the reactor cap using a number 8 wrench to unscrew the nuts of the height adjustable ports and screw them tight by hand. For exact locations of the ports see Fig. 3. The relative positions of the various ports are important in view of homogeneity, reproducibility and reliability of sensor readings. Make sure to use the tubes that are cut at an angle for inputs and for outputs the straight cut tubes. This is important because small additions of solution are more easily released from the tubes that are cut at an angle.
3. Now you can screw the cap on to the bottle and put a stirring rod in the reactor, don't forget to also add the centralizing disk between the reactor vessel and the screwcap. Then add demineralized water to the exact reactor volume with which you intend to run. Usually 60 ml is a handy volume to run with for the 100 ml reactor vessels, because it is easy volume to convert dilution rates. For example, if one runs at a flux-rate of 1 ml/per minute (with a working volume of 60 ml) this corresponds to a dilution rate of $1 \text{ v}\backslash\text{h}$.
4. Now the sensors should be calibrated (see section pH-sensor and O_2 -sensor).
5. Next you carefully insert the sensor(s). Crudely adjust their height so that they are deep enough that they sit into the liquid yet do not touch the stirring magnet on the bottom!
6. Then the bioreactor is placed on the magnet stirrer and the magnet stirrer is set at the speed of running (usually at 2/3 of the capacity). Now you will see that the volume of the water goes a bit down at the center and up at the sides. Also temporarily attach the gas input to the bioreactor so that gas is coming into the reactor and you can see that it works.
7. Next it is time to finely adjust the height of both sensors and the inputs and outputs. Sensors should be set at a depth of more or less 3 mm above the stirring rod when it is stirring. Now the efflux port can be adjusted. Unscrew the efflux port a little bit and adjust the height so that the tube just touches the water then screw the port tight (not too tight !). This will make sure that the volume of the reactor will always be at 60 ml. Now the bioreactor is ready for Autoclaving.

After connecting all of the tubing

1. The bioreactors should be autoclaved in the big autoclave on the 4th floor lab kitchen. The system is autoclaved for 20 min at 120°C. Make sure that enough water is added to the autoclave before starting the autoclave. An entire autoclaving session approximately takes 3-4 hours.
2. After removing the bioreactor frame from the autoclave add the stirring magnets into the frame sockets and put the frame into the preheated (30°C) waterbath. The tubing can then be placed back into the tubing guides and the peristaltic pumps, while the bioreactors further cool down. Do not yet secure the tubing cassettes into the pump.
3. For anaerobic running of the bioreactor the speed should be set at 1/3 of the way. (four o'clock position).

4. Connect the pH and possibly the DO sensor to the controllers. If a DO sensor is used it should be calibrated after autoclaving (see protocol).
5. Connect the gas inputs to the flowmeter. Check if the flowmeter indicates the correct gas flow 20 ml/min N₂ and additionally 100 ml/min air if the bioreactor is used for aerobic growth and adjust if necessary.
6. If you intent to use the pinch valves, check if the pinch valve power box is on.
7. Now connect the medium bottle, the efflux (waste) bottle, the acid/base and if necessary the migration tubing between bioreactors.
8. Now the tubing cassettes can be clicked into the pumps.
9. If necessary now perform the calibration of the DO sensor.
10. Start up the Labview application on the pc.
11. Activate the pump at high speed to wash out the demiwater from the bioreactors for four full dilutions.
12. Turn off the medium pump now and inoculate each bioreactor with 5 ml overnight culture via the sample port on the efflux tube.
13. Activate the pH control on the Labview application. Then grow the culture overnight in batch modus and sample the culture via the efflux sampling port.
14. If the culture is dense enough and all parameters (gas flow temperature, pH) are constant the experiment can start in continuous culture(bioreactor) modus.

Gluing stoppers after autoclaving

Because the welded-stopper-tubes are more than 2 times the price of the glued-stopper-tubes, it might sometimes be more economic to use glued stoppers and re-glue the stoppers after autoclaving using glue (*ERGO 5861*) and a primer (*ERGO 5150*). For gluing instructions see the protocol.

1. After the autoclaving the old glue should have turned from transparent to somewhat opaque and yellowish, this indicates that the old glue can be removed. All old glue should be removed carefully without damaging the tubing or the stopper and this area should be cleaned with a bit of ethanol and subsequently dried.
2. Then apply the primer (*Ergo 5150*) with the brush to both the tube and the stopper. This primer is based on an additive which is dissolved in an organic solvent. It serves to prepare non-polar and hard-to-glue polymers such as PharMed or silicone tubing. Let the primer dry before you start gluing.
3. Apply the adhesive (*Ergo 5861*) directly from the bottle on the stopper and glue the stopper at the exact location where it was before. This is to make sure that the tensioning of the tube will now deviate from the original specifications. Work fast because the glue dries fast! After 30 min the glue will be completely dry and the tubes are ready to use.

channel	NI 9208 I	NI 9208 II	NI 9403	NI 9870 I	NI 9870 II	NI 9870 III	USB hub
1	ph 1	ph 13	pinchvalve 1	julabo 1	julabo 5	reglo digital 3	reglo ICC 1
2	ph 2	ph 14	pinchvalve 2	julabo 2	julabo 6	ipcn-24 1	reglo ICC 2
3	ph 3	ph 15	pinchvalve 3	julabo 3	reglo digital 1	ipcn-24 2	reglo ICC 3
4	ph 4	ph 16	pinchvalve 4	julabo 4	reglo digital 2	ipcn-24 3	reglo ICC 4
5	ph 5	ph 17	pinchvalve 5				reglo ICC 5
6	ph 6	ph 18	pinchvalve 6				reglo ICC 6
7	ph 7	ph 19	pinchvalve 7				extrenal hard disk
8	ph 8	ph 20	pinchvalve 8				
9	shield	shield	shield				
10	shield	shield	shield				
11	ph 9	ph 21	pinchvalve 9				
12	ph 10	ph 22	pinchvalve 10				
13	ph 11	ph 23	pinchvalve 11				
14	ph 12	ph 24	pinchvalve 12				
15	DO 17		pinchvalve 13				
16	DO 18		pinchvalve 14				
17	DO 19		pinchvalve 15				
18	DO 20		pinchvalve 16				
19			RSVD				
20			pinchvalve 17				
21			pinchvalve 18				
22			pinchvalve 19				
23			pinchvalve 20				
24			pinchvalve 21				
25			pinchvalve 22				
26			pinchvalve 23				
27			pinchvalve 24				
28	shield	shield	shield				
29	shield	shield	shield				
30							
31							
32							
33							
34							
35							
36							
37							

Scheme 1.

**Chemostat
modules**

**Fuse box
and UPS**

**Chemostat
computer**




Fig. 1

air-input

N₂-input

pH-controllers

pH pump

**Medium
bottles**

**Migration
pump**




Fig. 2

Medium pump

Magnet stirrer

Waterbaths

**Chemostat
modules**




Fig. 3







Fig. 4








Fig. 8




Fig. 9




Fig. 10




Fig. 11




Fig. 12




Fig. 13



Fig. 14




Fig. 15





Fig. 16



Sensor screw connector

Sensor connector cap

Sensor solution cap

Fig. 17




Fig. 18





Fig. 19

N₂-distribution point

Air-distribution point



N₂-input

N₂-output

Air-input

Air-output

Fig. 20




Fig. 21

Air-output




Fig. 22




Fig. 23




Fig. 24




Fig. 25



Fig. 26





Fig. 27

Hold-power
tuning knob



On/off switch

Fig. 28

power input




Input control cRIO
(RS232)

Output control
pinch-valves

Fig. 29

Control cable



Dual channel
pinch-valve

Manual switch

Activation
indicator

Fig. 30



Fig. 31




Fig. 32




Fig. 33

pH-sensors and ICC-pumps Pinch-valves O₂-sensors

Control overview

Data logging displays




Fig. 34

Reglo digital pump

IPCN-24 pump

Waterbath




Fig. 35




Fig. 36




Fig. 37

Input elements editor

Input elements

Available program lines




Fig. 38




Fig. 39

Input and output parameters




Fig. 40