



1

2 **Supplemental figure 2: Platelets internalization by DCs enhances with increasing platelet:DC**  
 3 **ratios.** Freshly isolated platelets were labelled with PKH and incubated 1-4 hours with DCs.  
 4 Subsequently, DCs were harvested, fixed and stained with HLA-DR and CD61. Subsequently,  
 5 platelet internalization by DCs was quantified over time using imaging flow cytometry. Differences  
 6 in internalization were determined after 4 hours incubation using one-way ANOVA with Tukey post  
 7 testing. \*  $p < 0.05$ , \*\*  $p < 0.01$ , \*\*\*  $p < 0.001$ , \*\*\*\*  $p < 0.0001$ .

8